

# Extracting Nucleic Acids from UK NEQAS LI Samples



Annie Tapley

# What are UK NEQAS LI samples?

- Cultured cell lines (commercial or research)
- Blood donated from consenting patients/diagnostic waste
- Blood products (whole blood/buffy coat) NHSBT
- Bulking agents
- Samples are shipped as lyophilised material or fresh blood (usually the Chimerism programme only)

# Lyophilisation

- We freeze dry our samples a bit like ‘space food’ (please don’t eat our samples!)



- This process stabilises our samples enough to be posted worldwide to our participants - Saldanha *et al* (2007)

# Sample Stability

- If samples cannot be processed promptly upon receipt, they should be stored at 2-8 °C and remain lyophilised.
- This will ensure their stability until rehydration which revives the cells ready for laboratory processing but also revives RNases
- This is why we ask for sample processing to be performed immediately upon rehydration

# Sampling Handling

There are minimal changes required when handling UK NEQAS LI lyophilised samples:

1. They will need rehydrating in water (injection grade or better) - swirl and leave for a minute
2. There is no need for the red cell lysis step present in a lot of extraction methods (even if our samples are red in colour, lyophilisation will have lysed any red cells present)



# Sample Handling

3. Spin column style methods of DNA/RNA extraction often have a cell number limit
4. Sample cell numbers are provided in the accompanying paperwork - only approximate due to cell loss during lyophilisation/polyploidy/highly expressing
5. Dilute appropriately with PBS

# RNA extraction

- Acid guanidinium thiocyanate-phenol-chloroform extraction is the method used in-house by UK NEQAS LI to extract RNA
- It can be added directly to the lyophilised cells but will remain too viscous to do anything with unless they are left alone for at least five minutes
- Gentle but persistent pipetting will produce a homogenous solution ready for processing

# RNA extraction

- In-house extractions usually yield 40µl of RNA at concentrations of 600-1000ng/µl. More than enough template to synthesise cDNA from
- We have purchased a silica spin column based kit to get more hands on experience to aid in troubleshooting with participants



- UK NEQAS LI post cells instead of DNA/cDNA to assess **proficiency in extraction** as well as assay methodologies as this can have a fairly major impact on results
- We don't want this to prevent any participant from performing as well as you could so please don't hesitate to get in contact with any questions

# Don't Panic!

Everyone has 'off days'. No one is immune to human error/accidents

- Ask for repeat samples
- Ask for advice

We would be very grateful for your troubleshooting to be fed back to us so we can in turn help others



Questions?

[annie.tapley@ukneqasli.co.uk](mailto:annie.tapley@ukneqasli.co.uk)